




## EEG Notebooks: Mac installation instructions

### STEP 1: download the Miniconda installer for Mac

- Click the Python 3.6 file link on this website: <https://conda.io/miniconda.html>

## Miniconda

	 Windows	 MacOS	 Linux
Python 3.6	<a href="#">64-bit (exe installer)</a> <a href="#">32-bit (exe installer)</a>	<a href="#">64-bit (bash installer)</a>	<a href="#">64-bit (bash installer)</a> <a href="#">32-bit (bash installer)</a>
Python 2.7	<a href="#">64-bit (exe installer)</a> <a href="#">32-bit (exe installer)</a>	<a href="#">64-bit (bash installer)</a>	<a href="#">64-bit (bash installer)</a> <a href="#">32-bit (bash installer)</a>

- Press “command” and “spacebar” at the same time → type “Terminal” into the search bar → press the Terminal icon to open a new Terminal
- Enter the following commands in the Terminal:  

```
cd ~/Downloads/  
sh ./Miniconda3-latest-MacOSX-x86_64.sh
```
- You will be asked to type “yes” in the Terminal a couple times

## STEP 2: install required packages

- Run the following commands in the Terminal:

```
mv ~/eeg-notebooks ~/old-eeg-notebooks (if old version installed)
source ~/.bash_profile
cd ~
conda create -n nbmac python=3
source activate nbmac
conda install python=3.6
conda install git
git clone https://www.github.com/amandakeasson/eeg-notebooks
cd ~/eeg-notebooks
pip install -r requirements.txt
```
- For the “conda create” command, you will be asked to type “y” to confirm that yes, you wish to install the required packages
- The above packages may take ~5-10 minutes to install
- Ignore the warning that muselsl is not compatible with pygatt version 3.2.0
- Run the following commands in the Terminal:

```
mkdir ~/.jupyter/
jupyter notebook password
```
- You will then be asked to enter a password. Note that the password will not show up in the Terminal! After pressing <enter>, you will be asked to confirm your password
- Run the following commands in the Terminal:

```
python
import matplotlib
exit()
echo "backend: TkAgg" > ~/.matplotlib/matplotlibrc
```


## STEP 3: Working with eeg-notebooks

- Turn on your MUSE device
- Connect the Bluetooth dongle to the USB port of your Mac
- Run the following commands in the Terminal (in the nbmac environment)



```
source activate nbmac
cd ~/eeg-notebooks/notebooks
jupyter notebook
```
- Open **mac\_notebook.ipynb** in the browser
- The rest of the instructions for running experiments and analyzing data are in the mac\_notebook.ipynb notebook

## STEP 4: Shutting down eeg-notebooks

- Close the tab for the experiment notebook (e.g. mac\_n170.ipynb) in your browser
- In the jupyter notebook “Home” tab:
  - Click the checkbox beside the notebook(s) that is/are running (they will be green)
  - Click the orange “Shutdown” button that will appear after you click the checkbox
- In your terminal, press Ctrl + C, then press “y” and Enter when prompted to finish shutting down the jupyter notebook kernel

 jupyter Quit Logout

Files Running Clusters

Duplicate Shutdown View Edit  Upload New ▾ 

1 ▾ /

	Name ▾	Last Modified	File size
<input type="checkbox"/>	stimulus_presentation	a month ago	
<input type="checkbox"/>	utils	a month ago	
<input checked="" type="checkbox"/>	mac_N170.ipynb	Running a month ago	260 kB
<input type="checkbox"/>	N170.ipynb	a month ago	260 kB
<input type="checkbox"/>	P300.ipynb	a month ago	281 kB
<input type="checkbox"/>	Raw_EEG.ipynb	a month ago	1.42 kB
<input type="checkbox"/>	SSVEP.ipynb	a month ago	11.4 kB
<input type="checkbox"/>	_init_.py	a month ago	0 B
<input type="checkbox"/>	mac_n170.py	9 days ago	1.7 kB
<input type="checkbox"/>	run_experiment.py	a month ago	1.89 kB
<input type="checkbox"/>	test1.csv	a month ago	1.89 MB

```
(neurobrite_test2) ae: ~/code0ther/eeg-notebooks/notebooks $ jupyter notebook
[I 15:50:13.592 NotebookApp] Serving notebooks from local directory: /Users/amandae/code0ther/eeg-notebooks/notebooks
[I 15:50:13.593 NotebookApp] 0 active kernels
[I 15:50:13.593 NotebookApp] The Jupyter Notebook is running at:
[I 15:50:13.593 NotebookApp] http://localhost:8888/?token=6274037aa2845b5d5f34d5c5cb5c117c8353b5242cdba7fd
[I 15:50:13.593 NotebookApp] Use Control-C to stop this server and shut down all kernels (twice to skip confirmation).
[C 15:50:13.645 NotebookApp]

Copy/paste this URL into your browser when you connect for the first time,
to login with a token:
http://localhost:8888/?token=6274037aa2845b5d5f34d5c5cb5c117c8353b5242cdba7fd&token=6274037aa2845b5d5f34d5c5cb5c117c8353b5242cdba7fd
[I 15:50:17.356 NotebookApp] Accepting one-time-token-authenticated connection from ::1
[W 15:50:32.422 NotebookApp] Notebook mac_N170.ipynb is not trusted
[I 15:50:35.094 NotebookApp] Kernel started: cd633cec-206d-479e-b6a8-a8812c22452f
[W 15:50:45.598 NotebookApp] Timeout waiting for kernel_info reply from cd633cec-206d-479e-b6a8-a8812c22452f
[I 15:50:49.221 NotebookApp] Adapting to protocol v5.1 for kernel cd633cec-206d-479e-b6a8-a8812c22452f
[I 15:53:02.251 NotebookApp] Starting buffering for cd633cec-206d-479e-b6a8-a8812c22452f:e949634b178d457a83d8209c68505953
[I 15:55:11.334 NotebookApp] Kernel shutdown: cd633cec-206d-479e-b6a8-a8812c22452f

^C[I 15:56:32.650 NotebookApp] interrupted
Serving notebooks from local directory: /Users/amandae/code0ther/eeg-notebooks/notebooks
0 active kernels
The Jupyter Notebook is running at:
http://localhost:8888/?token=6274037aa2845b5d5f34d5c5cb5c117c8353b5242cdba7fd
Shutdown this notebook server (y/[n])? y
[C 15:56:34.298 NotebookApp] Shutdown confirmed
[I 15:56:34.365 NotebookApp] Shutting down 0 kernels
(neurobrite_test2) ae: ~/code0ther/eeg-notebooks/notebooks $
```